PROBE SYSTEM FOR LIFE



MINI PS4L DATA SHEET

m1-PS4L-2 (50 mm) Manual Probe System

The SemiProbe Mini-PS4L family of probe systems are built using the same design concepts as our patented Probe System for Life (PS4L) Architecture. The Mini-Series is designed to address the requirements of small sample testing using a small footprint.

The Mini-PS4L family of configured probe systems

addresses basic research and characterization applications that require probing on samples that are 100 mm and less in size. Our conventional family of probe systems - Lab Assistant and Probe System for Life (PS4L) are an overkill in size, capabilities, and cost for small sample testing applications.

The Mini-PS4L addresses the following applications:

- Basic Research and Device Characterization on small samples - die, partial wafers, PCB, substrates, and packaged parts
- Ideal for universities, government labs and corporations requiring simple testing capabilities
- Small Footprint easily fits on a bench, desk, cart or in a glovebox
- Cost configured systems start at \$6000
- Flexibility multiple system and accessory configurations are available
- Field Upgrade Path new modules and accessories can be added when required



m-PS4L-4 Kit Features & Benefits

- 50 mm Vacuum chuck
- Two complete MA-8005 manipulators with adjustable probe arm faceplates, coaxial probe arms, magnetic bases, and probe
- Fixed Platen fits up to four manipulators
- Fixed Stage
- Adjustable Optics (Z)
- Accessories manipulators, probe arms, probe tips, optics, CCTV System, Packaged Part Holders
- 100% Credit to PS4L Purchase entire mini purchase price credited to PS4L purchase

Major Applications/Markets Served

Device Characterization, MEMS, Optoelectronics, HF/Microwave, Failure Analysis, Photovoltaic, Research, Material & Life Science and more

The Mini-PS4L Series of probe systems are sold as a kitted configuration. Everything you need to start probing is included in the kit - prober, accessories, and installation manual. It is designed to be fully operational within 30 minutes after opening the box. This unique modular design enables customers to acquire test capabilities that precisely match their requirements and budget. More important, as the environment or test conditions change, the Mini-PS4L can be field-upgraded to add new capabilities or completely exchanged for a larger PS4L system at full credit. With this design philosophy, customers realize substantial time and cost savings over traditional probe systems because they do not need to invest in a new platform when wafer size, levels of automation or test requirements change.




SPECIFICATIONS

432 mm x 392 mm x 300 mm (17.0" x 15.4" x 12") (W,H,L) - with optics Base is 300 mm x 300 mm (12" x 12")
8 Kg (18 lbs.)
Fixed - chuck is mounted to a post that is bolted to the base of the system
Vacuum, round with concentric vacuum rings, or High Frequency (vacuum holes with calibration
substrate holders). Chucks will handle die, partial wafers, and wafers up to 75 mm. In addition,
substrates, PCBs, packaged parts, waffle packs, and more can be held in place
Other Chuck Options - 75 mm, 100 mm
Nickel plated steel with concentric vacuum rings (standard), other plating materials available
Planarity: +/- 8 µm
Aluminum with stainless steel top
Manipulator fixation - magnetic (standard), vacuum (optional with vacuum manifold kit) and bolt down (with holes
Fixed
No movement in X and Y, Adjustable in Z
Digital Camera or Compound Zoom Tube
Power: AC 110/220V AC 50-60 Hz 20A
Vacuum: 23 Hg or -0.8 bar

Note: Data, specifications and pictures vary depending on probe system configurations and accessories

PHYSICAL DIMENSIONS



Copyright 2011 SemiProbe, Inc. All rights reserved. Probe System for Life and SemiProbe are registered trademarks. Lab Assistant, IRIS Wafer Inspection and Pilot control software are trademarks of SemiProbe, Inc. All other trademarks are the property of their own respective owners.

276 E. Allen Street, Winooski, VT 05404 P: 802-860-7000 info@semiprobe.com www.semiprobe.com

